

COSMO*Pen*TM

The New Standard in
MICRONEEDLING TECHNOLOGY




A PRODUCT OF
 **COSMO**FRANCE

COSMOPen is a state-of-the-art **medical grade** high speed **Microneedling** system.


Get the same results as a fractional laser **without** pain, downtime, expense, and thermal damage.

COSMOPen can be used to treat **acne** and **surgical scars** - including **keloids** - as well as **wrinkles**, **stretch marks**, **pore size**, and signs of **aging**. Additionally, by vertically piercing the skin, **COSMOPen** helps naturally **induce collagen production** and **treat pigmentation**.

COSMOPen's new collagen creation:




DermaRollers



DermaRoller Flaws:

- Excessive & undesired trauma
- Cross infection risks
- Too large for small areas
- Lacks accuracy & consistency for speed & depth

COSMO*Pen*



COSMO*Pen* Advantages:

- Up and down movements prevent excessive trauma
- Disposable pinhead - No risk of cross infection
- The small tip helps treat hard to reach areas
- Very accurate & consistent speed & depth

VS

Cordless VS Coded

It's all about Convenience VS Efficiency

When designing the CosmoPen, we chose to make a **coded** device because of the following points:

Cordless:

Other

Pros:

- Convenience
- Mobility

Cons:

- Inconsistent power delivery*
- Power can run out during procedure
- Less powerful
- Wears out over time
- Heavier – Battery weight



Coded:

COSMOPen

Pros:

- Consistent power delivery*
- Can operate all day long
- Does not lose power over time – Longer life
- Lighter – No battery

Cons:


- Less mobility
- Less convenient

***Let's get technical:**

The amount of energy available is directly related to the state of charge of the battery. A battery that is 100% charged will significantly outperform a battery that is 60% charged. It is crucial for your device to be performing consistently at the selected speed during the entire microneedling procedure.

GET VISIBLE RESULTS!

COSMOPen uses an exclusive disposable 12-pin needle-head, which allows the user to control the depth of penetration (0.25 – 2.50 mm) and speed of the needles, up to 110 cycles/second, resulting in almost 1300 microchannels per second, which is 20% faster than most first generation devices.



Before treatment



**After 4 treatments with
CosmoPen & CBD (HealMD Topical)**



Before treatment



**After 3 treatments with
CosmoPen & CBD (HealMD Topical)**

Courtesy of **Kian Karimi MD, FACS**



The trusted CBD oil formulated by Kian Karimi MD, FACS

HealMD's Topical formula features 3,000mg of CBD isolate in an organic hempseed oil carrier. Those are the only two ingredients! It has been tested to ensure there are no chemicals, bacteria, yeast, or mold to be sure it is safe for topical use during Microneedling treatments.

Properties:

Anti-bacterial

Sebostatic

Anti-inflammatory

Indications:

To treat acne, helps with hyperpigmentation, acts as a powerful antioxidant, has anti-cancer components, and more.






Purchase supplies online:
www.cosmofrance.net



Videos & Information on:
www.cosmopen.net



📍 2121 SW 3rd Ave, Ste 500, Miami, FL 33129

☎ +1 (305) 538 0110 📠 +1 (305) 675 2642

💻 www.cosmofrance.net



A PRODUCT OF
COSMOFRANCE

CPCTLGREV5